



Instruction Manual

FastDNA[®] SPIN Kit for Soil

*Rapid Isolation of PCR - Ready Genomic
DNA from Soil Samples Using the FastPrep[®] System*

- › One Call**
- › One Source**
- › A World of
Biotechnology
Reagents**



Size:
50 preps

Storage:
Ambient temperature (15 – 30°C)

Catalog # 6560-200

Revision # 6560-200-07DEC

www.mpbio.com

MP Biomedicals • 29525 Fountain Parkway • Solon, OH 44139 • tel: 1.800.854.0530 • fax: 1.800.334.6999



Lyse BIG and FAST with Interchangeable Adapters!


QuickPrep™ Adapter

24 x 2 ml samples



TeenPrep™ Adapter

12 x 15ml samples



HighPrep™ Adapter

48 x 2ml samples



BigPrep™ Adapter

2 x 50 ml samples



The FastPrep-24 instrument is delivered with the QuickPrep™ Adapter

www.mpbio.com

FastDNA[®] SPIN Kit for Soil

*Rapid Isolation of PCR-Ready Genomic DNA
from Soil Samples Using the FastPrep[®] System*

50 Preps

Storage:
Ambient temperature (15 – 30°C)

Catalog # 6560-200

Revision # 6560-200-07DEC

www.mpbio.com



TABLE OF CONTENTS

1. Introduction to the FastDNA® SPIN Kit for Soil and the FastPrep® Instruments.....	5
2. Kit Components and User Supplied Materials	5
2.1 FastDNA® SPIN Kit for Soil Components.....	5
2.2 User Supplied Materials	6
3. Important Considerations Before Use	6
3.1 Preparation of SEWS-M Wash Solution	6
3.2 Sample Lysis with the FastPrep® Instrument	6
4. Safety Precautions	6
5. Protocol	7
6. Example Data : DNA Isolation from various Soil Samples and Gel Electrophoresis.....	8
7. Recommended Reference Format for Publication	9
8. References	9
9. Related Products.....	10
10. Product Use Limitation & Warranty	11

1. Introduction to the FastDNA® SPIN Kit for Soil and the FastPrep® Instruments

The FastDNA® SPIN Kit for Soil quickly and efficiently isolates PCR-ready genomic DNA directly from soil samples in less than 30 minutes. Designed for use with the FastPrep® Instruments from MP Biomedicals, plant and animal tissues, bacteria, algae, fungi spores and other members of a soil population are easily lysed within 40 seconds. These benchtop devices use a unique, optimized motion to homogenize samples by multidirectional, simultaneous impaction with lysing matrix particles. FastPrep® Instruments provide an extremely quick, efficient and highly reproducible homogenization that surpasses traditional extraction methods using enzymatic digestion, sonication, blending, douncing and vortexing. Samples are placed into 2.0 ml tubes containing Lysing Matrix E, a mixture of ceramic and silica particles designed to efficiently lyse all soil organisms including historically difficult sources such as eubacterial spores and endospores, gram positive bacteria, yeast, algae, nematodes and fungi. Homogenization in the FastPrep® Instrument with Lysing Matrix E takes place in the presence of MT Buffer and Sodium Phosphate Buffer, reagents carefully developed to protect and solubilize nucleic acids and proteins upon cell lysis. These reagents work together to allow extraction of genomic DNA with minimal RNA contamination.

Following lysis, samples are centrifuged to pellet soil, cell debris and lysing matrix. DNA is purified from the supernatant with a silica-based GENECLEAN® procedure using SPIN filters. Eluted DNA is ready for PCR, restriction digest, electrophoresis and any other desired application.

2. Kit Components and User Supplied Materials

2.1 FastDNA® SPIN Kit for Soil Components

Lysing Matrix E	50x 2.0 ml tubes
Sodium Phosphate Buffer	60 ml
MT Buffer	8 ml
PPS Solution	25 ml
Binding Matrix	66 ml
SPIN Modules	50 each
Catch Tubes	50 each
Concentrated SEWS-M	12 ml
DES	20ml
BBS Gel Loading Dye	200 µl
User manual	1 each
MSDS	1 each
Certificate of Analysis	1 each



FastDNA® Spin Kit for Soil

2.2 User Supplied Materials

FastPrep® Instrument (see Section 9)
Microcentrifuge that can freely spin 2.0 ml tubes
Microcentrifuge tubes (2.0 ml and 1.5 ml)
Clean 15 ml tubes for DNA binding
Rotator or low-speed vortex

3. Important Considerations Before Use

3.1 Preparation of SEWS-M Wash Solution


The FastDNA® SPIN Kit for Soil contains a bottle with 12 ml of Concentrated SEWS-M Wash Solution. Before using this solution, add 100 ml of 100% ethanol and mark on the bottle label the date ethanol was added. Ensure that the bottle is securely closed to prevent evaporation, mix and store at room temperature.

3.2 Sample Lysis with the FastPrep® Instrument

The fill volume in the lysing matrix tube after the addition of the Sodium Phosphate and MT Buffers to the sample should allow sufficient air space in the sample tube for efficient FastPrep® Instrument processing. MP Biomedicals recommends using 500 mg of starting material as long as there is between 250 – 500 µl of empty space in the tube. Sample loss or tube failure may result from overfilling the matrix tube. The matrix tube caps must be secure, but not over-tightened, to prevent sample leakage. If the sample is too large for processing in a single tube, divide the sample and process using multiple tubes. MP Biomedicals' Lysing Matrix particles and tubes have been rigorously tested and validated in the FastPrep® Instrument. The use of other products with the FastPrep® Instrument is not recommended and may result in sample loss or instrument failure. A single 40 second run at a speed setting of 6.0 in the FastPrep® Instrument is sufficient to lyse almost all samples. If the user experimentally determines that additional processing time is required, the sample should be incubated on ice in the Lysing Matrix E tube for at least 2 minutes between successive FastPrep® Instrument homogenizations to prevent overheating the sample and tube.

4. Safety Precautions

Binding Matrix contains components that, when in contact with human tissue, may cause irritation. Wear personal protective equipment to prevent contact with the skin or mucous membranes (gloves, lab coat, and eye protection). Consult the enclosed Material Safety Data Sheet for additional details.



5. Protocol

1. Add up to 500 mg of soil sample to a Lysing Matrix E tube.

NOTE: See section 3.2 for important guidelines.

2. Add 978 μ l Sodium Phosphate Buffer to sample in Lysing Matrix E tube.
3. Add 122 μ l MT Buffer.
4. Homogenize in the FastPrep[®] Instrument for 40 seconds at a speed setting of 6.0.
5. Centrifuge at 14,000 x g for 5-10 minutes to pellet debris.


NOTE: Extending centrifugation to 15 minutes can enhance elimination of excessive debris from large samples, or from cells with complex cell walls.

6. Transfer supernatant to a clean 2.0 ml microcentrifuge tube. Add 250 μ l PPS (Protein Precipitation Solution) and mix by shaking the tube by hand 10 times.
7. Centrifuge at 14,000 x g for 5 minutes to pellet precipitate. Transfer supernatant to a clean 15 ml tube. NOTE: While a 2.0 ml microcentrifuge tube may be used at this step, better mixing and DNA binding will occur in a larger tube.
8. Resuspend Binding Matrix suspension and add 1.0 ml to supernatant in 15 ml tube.
9. Place on rotator or invert by hand for 2 minutes to allow binding of DNA. Place tube in a rack for 3 minutes to allow settling of silica matrix.
10. Remove and discard 500 μ l of supernatant being careful to avoid settled Binding Matrix.
11. Resuspend Binding Matrix in the remaining amount of supernatant. Transfer approximately 600 μ l of the mixture to a SPIN[™] Filter and centrifuge at 14,000 x g for 1 minute. Empty the catch tube and add the remaining mixture to the SPIN[™] Filter and centrifuge as before. Empty the catch tube again.
12. Add 500 μ l prepared SEWS-M and gently resuspend the pellet using the force of the liquid from the pipet tip.
NOTE: Ensure that ethanol has been added to the Concentrated SEWS-M.
See section 3.1.

FastDNA® Spin Kit for Soil

13. Centrifuge at 14,000 x g for 1 minute. Empty the catch tube and replace.
14. Without any addition of liquid, centrifuge a second time at 14,000 x g for 2 minutes to “dry” the matrix of residual wash solution. Discard the catch tube and replace with a new, clean catch tube.
15. Air dry the SPIN™ Filter for 5 minutes at room temperature.
16. Gently resuspend Binding Matrix (above the SPIN filter) in 50-100 µl of DES (DNase/Pyrogen-Free Water).
NOTE: To avoid over-dilution of the purified DNA, use the smallest amount of DES required to resuspend Binding Matrix pellet.
NOTE: Yields may be increased by incubation for 5 minutes at 55°C in a heat block or water bath.
17. Centrifuge at 14,000 x g for 1 minute to bring eluted DNA into the clean catch tube. Discard the SPIN filter. DNA is now ready for PCR and other downstream applications. Store at -20°C for extended periods or 4°C until use.

6. Example Data: DNA Isolation from various Soil Samples and Gel Electrophoresis



DNA from various soil samples extracted with the FastDNA® SPIN Kit for Soil. 20% of the DNA isolated from 500mg soil was loaded on a 1.2% agarose gel (0.5X TAE). Soil was taken from:
Lane 1: tomato pot; Lane 2: sludge;
Lane 3: sandy soil; Lane 4: under pine tree;
Lane 5: under palm tree; Lane 6: green garden;
Lane 7: Nile Lilly pot; Lane 8: lawn grass;
Lane 9: citrus tree; Lane 10: avocado tree. DNA ranges from 4-20 kb.

7. Recommended Reference Format for Publications

DNA was isolated from (specific sample) using the FastDNA® SPIN Kit for Soil and the FastPrep® Instrument (MP Biomedicals, Santa Ana, CA).

8. References

Soil -

Jacob Bælum et al. (2006). Appl. Envir. Microbiol., Vol 72: 1476 - 1486.

Sediment -

Tracy J. Mincer et al. (2005). Appl. Envir. Microbiol., Vol 71: 7019 - 7028.

Snow samples -

Takahiro Segawa et al. (2005). Appl. Envir. Microbiol., Vol 71: 123 - 130.

Feces -

Alice Layton et al. (2006). Appl. Envir. Microbiol., Vol 72: 4214 - 4224

Sediments and Soil -

Francis C.A. et al. (2005). PNAS. Vol 102: 14683 - 14688

Soil -

Akira Ando et al. (2005). Appl. Envir. Microbiol., Vol. 71: 7075-7082

9. Related Products

Description	Size	Catalog #
FastPrep® 24 Instrument	100-230V	6002-500
FastPrep® FP100A Instrument	100V	6001-100
FastPrep® FP120A Instrument	120V	6001-120
FastPrep® FP220A Instrument	220V	6001-220
FastDNA® Kit	100 preps	6540-400
FastDNA® SPIN Kit	100 preps	6540-600
FastRNA® Pro Soil-Direct Kit	50 preps	6070-050
FastRNA® Pro Soil-Indirect Kit	50 Preps	6075-050
FastRNA® Pro Red Kit (Yeast & Fungus)	50 preps	6035-050
FastRNA® Pro Green Kit (Plant & Animal)	50 preps	6045-050
FastRNA® Pro Blue Kit (Bacteria)	50 preps	6025-050
FastProtein™ Blue Matrix	50 preps	6550-400
FastProtein™ Red Matrix	50 preps	6550-600
Lysing Matrix E	50 x 2 ml tubes	6914-050
Lysing Matrix E	100 x 2 ml tubes	6914-100
Lysing Matrix E	500 x 2 ml tubes	6914-500

10. Product Use Limitation & Warranty


Unless otherwise indicated, this product is for research use only. Purchase of MP Biomedicals products does not grant rights to reproduce, modify, or repackage the products or any derivative thereof to third parties. MP Biomedicals makes no warranty of any kind, expressed or implied, including merchantability or fitness for any particular purpose, except that the products sold will meet our specifications at the time of delivery. Buyer's exclusive remedy and the sole liability of MP Biomedicals hereunder shall be limited to, at our discretion, no replacement or compensation, product credits, refund of the purchase price of, or the replacement of materials that do not meet our specification. By acceptance of the product, Buyer indemnifies and holds MP Biomedicals harmless against, and assumes all liability for, the consequence of its use or misuse by the Buyer, its employees or others, including, but not limited to, the cost of handling. Said refund or replacement is conditioned on Buyer notifying MP Biomedicals, Inc. within thirty (30) days of receipt of product. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by the Buyer of all claims hereunder with respect to said material(s).

FastDNA[®], FastRNA[®], FastPrep[®] and BIO 101[®] Systems are registered trademarks of MP Biomedicals, LLC.



Ready-to-use Protocols For DNA, RNA And Protein Isolation From Any Sample

- Rapid and reproducible sample lysis and purification process
- No cross contamination with the closed lysing matrix tubes
- Increased yields of high quality DNA, RNA and proteins
- Integrity and size of DNA, RNA and proteins are retained
- Nucleic acids and proteins are ready-to-use in downstream application



FastDNA® Kit and FastDNA® Spin Kit

Cat N° 6540-400 - Cat N° 6540-600 respectively (100 preps)

- Lyse and isolate DNA in less than 30 minutes
- Plant, animal, yeast, fungal and microbial samples
- No hazardous organic reagents required
- SPIN filters streamline silica handling (FastDNA Spin Kit)

FastDNA® Spin Kit for Soil

Cat N° 6560-200 (100 preps)

- Lyse and isolate DNA in less than 30 minutes
- Variety of soil and environmental sample types
- No hazardous organic reagents required
- SPIN filters streamline silica handling

FastRNA® Pro Blue Kit

Cat N° 6025-050 (50 preps)

- For use with gram positive and gram negative bacteria
- Lyse up to 10^{10} cells per 2ml tube
- Lysis and isolation with single-phase organic solution in less than 90 minutes

FastRNA® Pro Red Kit

Cat N° 6035-050 (50 preps)

- For use with yeast cells and fungal tissue
- Lyse up to 10^{10} cells per 2ml tube
- Lysis and isolation with single-phase organic solution in less than 90 minutes

FastRNA® Pro Green Kit

Cat N° 6045-050 (50 preps)

- For use with all plant and animal samples
- Lyse 50-100 mg tissue per 2ml tube
- Lysis and isolation with single-phase organic solution in less than 90 minutes

FastRNA® Pro Soil-Direct Kit and FastRNA® Pro Soil-Indirect Kit

Cat N° 6070-050 - Cat N° 6075-050 respectively (50 preps)

- Isolate RNA from soil samples (direct kit) and washed soil (indirect kit) in less than 2 hours
- Variety of soil and environmental sample types
- RNA protected during and after processing
- Humic acids reduced to allow uninhibited RT-PCR
- Includes additional reagents for even further purification if necessary
- SPIN filters streamline silica handling

FastProtein™ Blue Matrix

Cat N° 6550-400 (50 preps) - Cat N° 6550-500 (100 preps)


- Release of proteins from gram positive and gram negative bacteria in 40 seconds
- Protein extracts are ready for immediate electrophoresis or purification
- Ideal for optimizing induction conditions

FastProtein™ Red Matrix

Cat N° 6550-600 (50 preps) - Cat N° 6550-700 (100 preps)

- Release of proteins from yeast cells and fungi in 40 seconds
- Protein extracts are ready for immediate electrophoresis or purification
- Ideal for optimizing induction conditions

MP



TM

MP Biomedicals

Our family of companies includes:



Lyse BIG and FAST with Interchangeable Adapters!


QuickPrep™ Adapter

24 x 2 ml samples




TeenPrep™ Adapter

12 x 15ml samples



HighPrep™ Adapter

48 x 2ml samples



BigPrep™ Adapter

2 x 50 ml samples



The FastPrep-24 instrument is delivered with the QuickPrep™ Adapter

www.mpbio.com



Instruction Manual

FastDNA[®] SPIN Kit for Soil

Protocol Revision # 6560-200-07DEC

Worldwide Ordering and Technical Support

United States of America

Worldwide Headquarters

Tel: +1.440.337.1200

Toll Free Tel: 800.854.0530

Fax: +1.440.337.1180

Toll Free Fax: 800.334.6999

Europe

Toll Free Phone:

00800.7777.9999

Toll Free Fax: 00800.6666.8888

Australia

MP Biomedicals Australasia Pty Ltd

Tel: +61.2.9838.7422

Fax: +61.2.9838.7390

Belgium

MP Biomedicals

Tel: 02 466 00 00

Fax: 02 466 26 42

Canada

MP Biomedicals Canada

Tel: 888.362.5487

Fax: 514.935.7541

France

MP Biomedicals France

Tel: 03 88 67 54 25

Fax: 03 88 67 19 45

Germany

MP Biomedicals

Phone: 0800 426 67337

Fax: 0800 629 67337

Japan

MP Bio Japan K.K.

Tel: 03-3808-2102

Toll Free Tel: 0120.788.020

Fax: 03-3808-2401

The Netherlands

MP Biomedicals Netherlands

Tel: 0800-0227416

Fax: 0800-0227489

Poland

MP Biomedicals Poland

Tel: +48.22.659.58.95

Fax: +48.22.658.45.05

Serbia

MP Global d.o.o.

Tel: +381.11.2622.945

Fax: +381.11.2623.373

Singapore

MP Biomedicals Singapore

Tel: 65.6775.0008

Fax: 65.6775.4536

Switzerland

MP Biomedicals Switzerland

Tel: 061 271 0007

Fax: 061 271 0084

United Kingdom

MP Biomedicals UK

Tel: 0800 282 474

Fax: 0800 614 735

www.mpbio.com

MP Biomedicals • 29525 Fountain Parkway • Solon, OH 44139 • tel: 1.800.854.0530 • fax: 1.800.334.6999

